




Endoscopy Ultra Sound Aspiration needle / EUS

is used for the sampling of the submucosal
gastrointestinal space and lesions




Endoscopy Ultra Sound aspiration needle



EUS needle

- optimal view of the cannula due to grooved structure underneath the needle tip
- 5 cm adjustable tube position, suitable for all biopsy channels
- kink resistant mandrin rod
- 8 cm adjustable, graduated needle pusher
- smooth-running, kink resistant plastic insertion tube
- highly flexible, kink resistant stainless steel cannula with facet grinding
- tube length adjustable



Article numbers **EUS needles**

Article-No.		Diameter mm / Fr	Length cm	Needle Diameter Gauge	PU
90TU 19-A	Mandrin with ball tip	1,75 / 5,2	142	19	2
90TU 22-A	Mandrin with ball tip	1,75 / 5,2	142	22	2
90TU 25-A	Mandrin with ball tip	1,75 / 5,2	142	25	2

We reserve the right to make technical changes and to deviate from the depicted images.

