


Cleaning Brush for Endoscope

Sterile cleaning brush for sampling the working channels of endoscopes

- Double-sided, sterile
- One-sided, sterile



STERILE EO

MTW[®]
Endoskopie Manufaktur



MTW-Endoskopie W. Haag KG
Goldsbergstr. 18, 46487 Wesel, Germany
Tel: +49 2803 9120-0
info@mtw-endoskopie.com
www.mtw-endoskopie.com

MD Cleaning Brush for Endoscope


Model: GM-CB-A-0505-18-2300

REF CB-A-0505-18-2300

LOT CB25111302

EN 2025-11-13

EXP 2028-11-12



de | Reinigungsbürste für Endoskop
fr | Brosse de Nettoyage pour Endoscope
es | Cepillo de Limpieza para Endoscopio
it | Spazzola di Pulizia per Endoscopio
nl | Reinigingsborstel Voor Endoscoop
no | Rengjøringsbørste for Endoskop
pt | Escova de Limpeza Para Endoscópio



(01) 0 6971890 56477 1



(17) 281112 (10) CB25111302

Cleaning Brush for Endoscope / Reinigungsbürste für Endoskop

MTW[®]
Endoskopie Manufaktur

REF CB-A-0505-18-2300 **EXP** 2028-11-12 **LOT** CB25111302

UDI (01) 0 6971890 56477 1 (17) 281112 (10) CB25111302

Cleaning Brush for Endoscope / Reinigungsbürste für Endoskop

MTW[®]
Endoskopie Manufaktur

REF CB-A-0505-18-2300 **EXP** 2028-11-12 **LOT** CB25111302

UDI (01) 0 6971890 56477 1 (17) 281112 (10) CB25111302

Cleaning Brush for Endoscope / Reinigungsbürste für Endoskop

MTW[®]
Endoskopie Manufaktur

REF CB-A-0505-18-2300 **EXP** 2028-11-12 **LOT** CB25111302

UDI (01) 0 6971890 56477 1 (17) 281112 (10) CB25111302


Cleaning Brush for Endoscope / Reinigungsbürste für Endoskop

MTW[®]
Endoskopie Manufaktur

REF CB-A-0505-18-2300 **EXP** 2028-11-12 **LOT** CB25111302

UDI (01) 0 6971890 56477 1 (17) 281112 (10) CB25111302

STERILE EO



Order no.	Article no.	Description	Brush diameter mm	Working channel mm	Diameter mm	Length cm	PU
CB-A-0303-18-2300 x20	CB-A-0303-18-2300	double-sided	3,0	2,0-2,6	1,8	230	20
CB-A-0505-18-2300 x20	CB-A-0505-18-2300	double-sided	5,0	2,7-4,6	1,8	230	20
CB-D-03-18-2300 x20	CB-D-03-18-2300	one-sided	3,0	2,0-2,6	1,8	230	20
CB-D-05-18-2300 x20	CB-D-05-18-2300	one-sided	5,0	2,7-4,6	1,8	230	20

