



Papillotome

„Safe-Cut“, partially insulated cutting wire





Proximal 3-lumen connector

A 3-lumen connector allows for the simultaneous administration of contrast medium while the guidewire is in place.



Papillotome

„Safe-Cut“, partially insulated cutting wire

Order no.	Article no.	Special features	Feather strengthened	For guide wire inch	Cutting length mm	Diameter mm	Length cm	PU
99028141 x5	99028141	8 mm guiding tip	yes	0.035	30, partially insulated	2,3 > 1,8	215	5

We reserve the right to make technical changes and to deviate from the depicted images.

