

Disposable Ultrasound Biopsy Needles

EUS-FNB-Biopsy Needle (Stainless steel)

EUS-FNA-Biopsy Needle (Stainless steel)

EBUS-TBNB-Biopsy Needle (Stainless steel)

EBUS-TBNA-Biopsy Needle (Stainless steel)

EUS-FNB-Biopsy Needle (Nitinol)

EUS-FNA-Biopsy Needle (Nitinol)



FNB


FNA

Options:

- Highly flexible 100% Nitinol needles
- Stainless steel

The nitinol needle core wire





detailed information

Disposable Ultrasound Biopsy Needles

EUS-FNB-Biopsy Needle (Stainless steel)

EUS-FNA-Biopsy Needle (Stainless steel)

EBUS-TBNB-Biopsy Needle (Stainless steel)

EBUS-TBNA-Biopsy Needle (Stainless steel)

EUS-FNB-Biopsy Needle (Nitinol)

EUS-FNA-Biopsy Needle (Nitinol)

- **Optimal visibility of the cannula** thanks to the “hole structure” below the needle tip
- Kink-resistant Nitinol mandrel
- Smooth-running, kink-resistant plastic insertion tube
- Highly flexible, kink-resistant variants: stainless steel & **100% Nitinol** needles

		Articel no.	Description			Needle Ø Gauge	Need- le mm	Ø mm	Length mm	PU
OLYMPUS FUJIFILM PENTAX		UN-A-19-11-80-16-1375	FNA Stainless steel EUS	FNA		19	1,1	1,6	1375	1
		UN-A-22-07-80-16-1375	FNA Stainless steel EUS	FNA		22	0,7	1,6	1375	1
		UN-B-22-07-80-16-1375	FNA Stainless steel EUS	FNB		22	0,7	1,6	1375	1
		UN-B-22-07-80-22-1375	FNA Stainless steel EUS	FNB		22	0,7	2,2	1375	1
		UN-D-22-07-80-16-1375	FNB Nitinol EUS	FNB		22	0,7	1,6	1375	1
		UN-C-22-07-80-16-1375	FNA Nitinol EUS	FNA		22	0,7	1,6	1375	1
		UN-D-22-07-80-22-1375	FNB Nitinol EUS	FNB		22	0,7	2,2	1375	1
OLYMPUS FUJIFILM		UN-A-19-11-50-16-715	EBUS-TBNA, Stainless steel	FNA		19	1,1	1,6	715	1
		UN-A-22-07-50-16-715	EBUS-TBNA, Stainless steel	FNA		22	0,7	1,6	715	1
		UN-B-22-07-50-16-715	EBUS-TBNB, Stainless steel	FNB		22	0,7	1,6	715	1
		UN-C-19-11-50-16-715	EBUS-TBNA, Nitinol	FNA		19	1,1	1,6	715	1
		UN-C-22-07-50-16-715	EBUS-TBNA, Nitinol	FNA		22	0,7	1,6	715	1
		UN-D-22-07-50-16-715	EBUS-TBNB, Nitinol	FNB		22	0,7	1,6	715	1
PENTAX		UN-A-22-07-50-16-670	EBUS-TBNA, Stainless steel	FNA		22	0,7	1,6	670	1
		UN-A-19-11-50-16-670	EBUS-TBNA, Stainless steel	FNA		19	1,1	1,6	670	1
		UN-B-22-07-50-16-670	EBUS-TBNB, Stainless steel	FNB		22	0,7	1,6	670	1

We reserve the right to make technical changes and to deviate from the depicted images.

